

# STERILABEL label/barcode printer



SERIAL NUMBER: 16019 - CYCLE No.: 037  
 134 DENTAL PIECES S1  
 14/09/09 15:03 - 15:16  
 USER Factory - Cycle completed  
 USE BY: 14/10/09

## STERILABEL

	M.U.	
Labels	n.	300 labels/roll
Printing data		serial no - cycle no - cycle type - date - starting/finishing time - user - cycle result - packaging expiry date (if setted)
Interface		serial RS232
Supply	VDC/VA	9-48 VDC / 11 VA (external supplier)
Dimensions	mm	146 x 88 x 65
Weight	g	350 g (with labels roll)

Technical data are subject to change without notice

STERILABEL is the thermal label printer which can record the sterilization process in order to **recognize** materials and check the process **traceability**: it can be really easily used. Up to 100 simple and **precut labels** can be print for every single sterilization cycle. Also **BARCODE** can be print.