

# Completely equipped

## Supplement set for VistaScan

Everything you could want in materials for using VistaScan models is provided in our supplement sets.



### Supplement Set VistaScan Plus / Perio Plus

Image plate Size 1, Size 3, 2 pieces each  
Image plate Size 4, 1 piece  
Foil cassette Size 1, Size 3, Size 4, 1 piece each  
Light protection cover Size 1, Size 3, Size 4, 100 pieces each  
Bit protection Size 4, 100 pieces  
**2130-007-00 VistaScan Plus**  
**2140-007-00 VistaScan Perio Plus**

### Supplement Set VistaScan Mini Plus

Image plate Size 1, Size 3, 2 pieces each  
Image plate Size 4, 1 piece  
Light protection cover Size 1, Size 3, Size 4, 100 pieces each  
Bite protection Size 4, 100 pieces  
**2141-007-00**

# Highest image resolution

## Image Plate Plus

Thanks to an image resolution of 22 lp/mm, Image Plate Plus is superior X-ray film in recognising details – and therefore superior to that of the competition. The reusable image plates are especially thin and flexible and therefore comfortable for the patients.



## Our intraoral image plate offer

### Image Plate Plus Size 0

Format 2 x 3 cm, 2 pieces  
**2130-040-50**  
**2130-040-60 with copper point and ID**

### Image Plate Plus Size 1

Format 2 x 4 cm, 2 pieces  
**2130-041-50**

### Image Plate Plus Size 2

Format 3 x 4 cm, 4 pieces  
**2130-042-50**  
**2130-042-60 with copper point and ID**

### Image Plate Plus Size 3

Format 2.7 x 5.4 cm, 2 pieces  
**2130-043-50**

### Image Plate Plus Size 4

Format 5.7 x 7.5 cm, 1 piece  
**2130-044-50**

P007-597-02/HK-Dom/3\_01/07Y Subject to change without notice

DÜRR DENTAL AG  
Höpfigheimer Strasse 17  
74321 Bietigheim-Bissingen  
Germany  
www.duerr.de  
info@duerr.de

# VistaScan – image plate accessories



COMPRESSED AIR  
SUCTION  
IMAGING  
DENTAL CARE  
HYGIENE



# Safe application

## Optimal protection

The light protection cover for image plates for protecting against exposure-to light and contamination are made of health-tested material, are 100% hygienically sealed and have thin-edges and a recess on each side for left and right handed people. Studies have proven the high light-proof capabilities.



### Light protection cover Size 0

Protective sleeves for image plate format 2 x 3 cm, 100 pieces  
**2130-070-00/2130-070-50\***

### Light protection cover Size 1

Protective sleeves for image plate format 2 x 4 cm, 100 pieces  
**2130-071-00**

### Light protection cover Size 2

Protective sleeves for image plate format 3 x 4 cm, 300 pieces  
**2130-072-00/2130-072-50\***

### Light protection cover Size 3

Protective sleeves for image plate format 2.7 x 5.4 cm, 100 pieces  
**2130-073-00**

### Light protection cover Size 4

Protective sleeves for image plate format 5.7 x 7.5 cm, 100 pieces  
**2130-074-00 VistaScan Plus/Mini Plus**  
**2140-074-00 VistaScan Perio Plus**

### Bite protection Size 4

Bite protection for image plate Size 4 (occlusal), 100 pieces  
**2130-074-03**

# Easy handling



## Image plate holder system

The image plate holders by Dürr Dental offer excellent protection for the image plate thanks to their de-burred and rounded edges. There are holders for all applications to ensure precise X-rays images.

### Image plates and film holder system set

Sight-ring anterior and bite wing, 2 pieces  
Sight-ring posterior and endo, 2 pieces  
Bite plates anterior, 6 pieces  
Bite plates posterior, 6 pieces  
Bite plates Bite wings 3 x 4 cm, 2 pieces  
Bite plates Bite wings 2 x 3 cm, 2 pieces  
Guide rod anterior, 1 piece  
Guide rod posterior, 2 pieces  
Guide rod Bite wing, 1 piece

**2130-981-50**

### Image plates and film holding system retrofitted for endo recording

Guide rod endo, 1 piece  
Bite plate endo Left anterior, 1 piece  
Bite plate endo Left posterior, 1 piece  
Bite plate endo Right anterior, 1 piece  
Bite plate endo Right posterior, 1 piece

**2130-981-51**

### Bite plate Bite wing

2.7 x 5.4 cm, 8 piece

**2130-981-17**

# Disinfection & Cleaning

## FD 350 disinfectant wipes

Especially suitable for wipe disinfection and cleaning of medical products. For example sleeves of intraoral foils and films, hand-pieces and angle-pieces etc.

### FD 350 disinfection wipes

12 x container of 90 wipes each

**CDF350C0140**

### Refill pack

2 x 6 refill packs of 90 wipes

**CDF350C1240**



## Image plate cleaning cloth

The VistaScan image plate cleaning wipes clean the film intensively and also disinfects. The benefit is a lasting, optimal image quality.

### Image plate cleaning cloth

10 wipes

**CCB351B1001**



\*also available in white